



PURE PLANET PAPER CUPS

HIGHLY QUALITATIVE DENTAL CUP



CONFIGURATION OVERVIEW


- Material:** Water based coated Paper
- General characteristics:**
- The product meets the requirements for food contact
 - Compostable, made of eco-friendly virgin wood pulp
 - Produced according to the following standards = ISO14001:2015, ISO 45001:2018, ISO 9001:2015.
 - Water based compostable film 5-10 g/m2 (NO PE)

Colors: White star

TECHNICAL SPECIFICATIONS

Weight	210 gsm paper
Volume	200 ml = L01
Measures	73x50x78 mm (top diameter/bottom diameter/height)

PRODUCT VARIATIONS

Color code	WH
Color	White star
	

PACKAGING

Pcs/carton:	20x50 pcs
Generic Art-No.	CUPAY50-AKZ-WH-L01
Carton size:	37x30.5x31 cm
Carton weight:	4,6 Kg
Carton/layers:	6
Layer/pallet:	6
Carton/pallet:	36
Pallet size:	80x120x206
Pallet weight:	180.6 Kg

STORAGE / WASTAGE

Store in a dry, dark and airy space away from heat sources.

!ATTENTION! It is recommended to use this product only with water. Adding chemicals could damage the structure of the product. If used with other substances other than water it is suggested to use the solution instantly and to discard the cup immediately after use. Akzenta won't be held responsible for any damages resulting from an improper use of the product.

WARNING : This product may contain latex.