

# TOP TOUCH PLUS NITRILE ULTRASOFT GLOVES

## SINGLE USE DISPOSABLE NITRILE EXAMINATION GLOVES

### CONFIGURATION OVERVIEW

<b>Features</b>	Cuff strengthened by a rolled cuff Easy donning and removal Excellent tactile sensibility (Micro textured fingertips) Suitable for food handling purposes
<b>Shape:</b>	Ambidextrous
<b>Material:</b>	Nitrile
<b>Type:</b>	Non-Sterile Powder-free Latex-free

### TECHNICAL SPECIFICATIONS

#### Physical Properties:

Elongation at Break (%)	≥500%
Force Break (N)	≥6N

Location of Thickness Measurement	Single Wall (mm)
Finger Tip	≥ 0.08
Palm	≥ 0.05
Cuff	≥ 0.04

<b>Available Length:</b>	≥ 240 mm	<b>Width :</b> 70 ( XS), 80 (S), 95 (M), 110 (L)
<b>Weight M size:</b>	3.5 g ± 0.2 g	
<b>Shelf life:</b>	3 years	

### QUALITY ASSURANCE

#### Quality certification:

Class I Medical Device (under classification MDR, Annex VIII, rule 5)  
Conform to EU standard EN 455 part 1:2020, part 2:2015, part 3:2015, part 4:2009  
Conform to EN 1041:2008+A1:2013, EN ISO 15223-1:2017, EN ISO 14971:2012, EN 16523-1:2015  
Produced under ISO 9001:2015 and 13485 :2016  
Conform to ISO 10993-10:2010, 10993-5:2009 and 10993-1:2009  
Standard and technical specifications applied : EN420:2003+A1:2009  
Conform to EU REG. 2016/425 standard EN 374 PART 1:2016, part 2:2014,part 4:2013 and part 5:2016

### PRODUCT VARIATIONS

**Available Sizes:** XS S M L

#### PACKAGING

<b>Pcs/Dispenser:</b>	100
<b>Dispenser/Carton:</b>	10
<b>Carton Size:</b>	31x22.5x21
<b>Carton weight:</b>	XS: 3.5 Kg. S: 3.75 Kg. M: 4.25 Kg. L: 4.75 Kg.
<b>Pallet size:</b>	cm 80x120x209
<b>Pallet weight :</b>	XS : 330 Kg. S : 352.5 Kg. M : 397.5 Kg. L : 442.5 Kg.
<b>Carton / layer:</b>	10
<b>Layer / pallet:</b>	9
<b>Carton / pallet:</b>	90

Art- No.	GLOBN51-AKZ- LA	GLOBN51-AKZ- WH
Color	Lavender	White
		

