



STYLE NITRILE BLACK GLOVES

SINGLE USE DISPOSABLE NITRILE EXAMINATION GLOVES



CONFIGURATION OVERVIEW

- Features**
- Cuff strengthened by a rolled cuff
 - Easy donning and removal
 - Excellent tactile sensibility (Micro textured fingertips)
 - Suitable for food handling purposes
- Shape:** Ambidextrous
- Material:** Nitrile
- Type:** Non-Sterile
Powder-free
Latex-free

TECHNICAL SPECIFICATIONS

Physical Properties:

Elongation at Break (%)	≥500%
Force Break (N)	≥6N
Location of Thickness Measurement	Single Wall (mm)
Finger Tip	≥ 0.08
Palm	≥ 0.05
Cuff	≥ 0.04

- Available Length:** ≥ 240 mm **Width :** 70 (XS), 80 (S), 95 (M), 110 (L)
- Weight M size:** 3.5 g ± 0.2 g
- Shelf life:** 3 years

QUALITY ASSURANCE

Quality certification: Class I Medical Device (under classification MDR, Annex VIII, rule 5)
 Conform to EU standard EN 455 part 1:2020, part 2:2015, part 3:2015, part 4:2009
 Conform to EN 1041:2008+A1:2013, EN ISO 15223-1:2021, EN ISO 14971:2019, EN 16523-1:2015
 Produced under ISO 9001:2015 and 13485 :2016
 Conform to ISO 10993-10:2021, 10993-5:2009 and 10993-1:2018
 Standard and technical specifications applied : EN420:2003+A1:2009
 Conform to EU REG. 2016/425 standard EN 374 PART 1:2016, part 2:2014,part 4:2013 and part 5:2016

PRODUCT VARIATIONS

Available Sizes: XS S M L

Art-No.	GLOBN51-AKZ-BK
Color	Black

PACKAGING

- Pcs/Dispenser:** 100
- Dispenser/Carton:** 10
- Carton Size:** 31x22.5x21 cm
- Pallet size:** cm 80x120x209
- Carton weight:** XS : 3.50Kg. S: 3.75Kg. M: 4.25 Kg. L : 4.75Kg
- Pallet weight :** XS : 330 Kg. S : 352.5 Kg. M : 397.5 Kg. L : 442.5 Kg.
- Carton / layer:** 10
- Layer / pallet:** 9
- Carton / pallet:** 90



Akzenta International SA
Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
fax. +41 (0) 91 682 85 20

www.akzenta.com
e-mail: business@akzenta.com