



MICRO APPLICATORS

DISPOSABLE MICRO APPLICATORS



CONFIGURATION OVERVIEW

- Material:** Medical level Polypropylene handle, Nylon Brush
Type: Non-Sterile
Description:
- Ideal for application of sealants, bonding agents, etchants, hemostatic solutions, cavity liners and disclosing solutions.
 - Non-absorbent, non-linting fibers suspend solutions as small as 1/8 of a drop without spilling, dripping and waste.
 - Bendable neck "stays put" to facilitate the placement of material into difficult areas.
 - Flexible plastic handle for pinpoint accuracy application of any liquids
 - Class 1 Medical Device according to EU MDR 20017/745 according to rule 5 Annex VIII
 - Harmonised to : EN ISO 14971:2012, EN ISO 15223-1:2016, ISO 10993-1:2018, EN ISO 10993-5:2009, EN ISO 10993-10:2013, EN 1041:2008+A1:2013

THECNICAL SPECIFICATION

| Size | Regular | Fine | Superfine |
|---|---------|--------|-----------|
| Length of the applicator (including the fiber tip) (mm) | 102 | 101.45 | 100.2 |
| Diameter of the fiber tip (mm) | 3 | 1.9 | 1.5 |
| Length of the fiber tip (mm) | 3.5 | 3.5 | 2.6 |

PRODUCT VARIATIONS

| Regular | Art-No. | MICAP25-AKZ-BL-REG | MICAP25-AKZ-LI-REG | MICAP25-AKZ-FG-REG | MICAP25-AKZ-OR-REG | MICAP25-AKZ-FU-REG | MICAP25-AKZ-PI-REG |
|---------|---------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|
| | Color | Med Blue | Lilac | Fresh Green | Orange | Fuchsia | Pink |

| Fine | Art-No. | MICAP25-AKZ.YE.FNE | MICAP25-AKZ-MI-FNE | MICAP25-AKZ-MC-FNE | MICAP25-AKZ-LU-FNE | MICAP25-AKZ-CE-FNE |
|------|---------|--------------------|--------------------|--------------------|--------------------|--------------------|
| | Color | Yellow | Mint Green | Magic Blue | Laguna Blue | Crispy cedro |

| Superfine | Art-No. | MICAP25-AKZ-WH-SFN |
|-----------|---------|--------------------|
| | Color | White |

PACKAGING

| | |
|---------------|-----------------|
| Pcs/dispenser | 100 |
| Dispenser/box | 4 |
| Boxes/carton | 40 |
| Carton size | 52x42x24.5 cm |
| Carton weight | 16 Kg |
| Carton/pallet | 33 |
| Pallet size | 80x120x235,5 cm |
| Pallet weight | 548 Kg |

SHELF LIFE : 5 years

STORAGE AND HANDLING

In dry, dark, dust-free and airy space, free of disturbing smells.

USAGE METHOD : Choose the suitable size of applicators to be coated with medicine to the affected area.



Akzenta International SA
 Via Giuseppe Motta, 24
 CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
 fax. +41 (0) 91 682 85 20

www.akzenta.com
 e-mail: business@akzenta.com