



TOP FLEX PLUS SALIVA EJECTOR



DISPOSABLE SALIVA EJECTOR

CONFIGURATION OVERVIEW

Material: Cannula: non-toxic polyvinyl chloride PVC with a zinc- or copper-coated iron wire insert
Smooth surface
Filtering part: Polyethylene (LDPE), while in the model with fixed hood is made from polyvinyl chloride (PVC) – removable
Filter with No. 6 openings of approx. 1 mm

Latex: Absent

Dangerous components classified by EC: Absent

Colors: White / Transparent / Yellow / Pink / Fuchsia / Red / Lilac / Burgundy / Fresh green/ Med blue / Dark green / Dark blue / Light blue / Orange / Mint green

TECHNICAL SPECIFICATIONS

Cannula

External diameter	6.5-6.6 mm
Thickness	1.0 + 1.2 mm
Length	145 mm or 125 mm
Metal wire	Same length as the cannula

Filtering Part

Length	14 mm
Outside diameter	9 mm

QUALITY ASSURANCE

Quality certification: ISO 9001:2015 and ISO 13485:2016, UNI CEI ISO 14971:2020, UNI CEI EN 15986:2011, UNI CEI ISO 15223-1:2017, skin sensitization test and irritation test ISO 10993-10:2010 and ISO 10993-1:2018, in vitro cytotoxicity ISO 10993-5:2009 ; UNI EN 1186-3:2003; ISO 15223-1:2017, CEI UNI EN 1041:2014

Conformity to guidelines: medical device of class IIA EU REG. 2017/745 Annex I

SHELF LIFE 5 Years

Article Number: SALEA03-AKZ...C44 (145mm) / SALBA03-AKZ...C45 (125mm)

Art-No.	WH	TR	YE	PI	FU	RE	LI	BU
Color	White	Transparent	Yellow	Pink	Fuchsia	Red	Lilac	Burgundy
Art-No.	FG	MB	DG	MC	LB	HO	EG	
Color	Fresh Green	Med blue	Dark green	Dark blue	Light blue	Orange	Mint green	

PACKAGING

Format	125mm/145 mm	Carton weight	4.6 Kg
Pcs/pack	100		
Pack/carton	10		
Carton size	29x18x20.5 cm		
Cartons/layer	16		
Layers/pallet	9		
Cartons/pallet	144		
Palette Size	80 x 120 x 200		



Akzenta International SA
Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 91 210 46 00
fax. +41 91 682 85 20

www.akzenta.com
e-mail: business@akzenta.com