



MERCATOR dermagel coated

Glove Type:	Natural rubber latex powder-free, examination and protective, non-sterile
Color:	Creamy
Shape:	Ambidextrous, gloves fitting either hand
Powder:	Powder-free
Protein content:	<50 µg/g,
Outer surface:	textured
Inner surface:	polymer coated
Cuff:	beaded
Min. length [mm]:	240
Thickness [mm] (finger /palm/cuff):	0,12 / 0,09 / 0,06
Force at break [N] (before /after ageing):	6N / 6N
AQL:	1.5
Classification:	Medical Device: CE class I (Regulation (UE) 2017/745) Personal Protective Equipment: CE2777 Category III (Regulation (UE) 2016/425) Type C (EN ISO 374-1)
Compliance with:	EN 455 (1-4), EN ISO 15223-1, EN 20417, EN ISO 374-1 (Typ C), EN ISO 374-2, EN ISO 374-4, EN ISO 374-5, EN 16523-1, EN ISO 21420, ISO 16604
Pictograms:	
Sizes:	 Size Color Coding
REF No:	XS - RD10006001 S - RD10006002 M - RD10006003 L - RD10006004 XL - RD10006005
Logistic information:	Quantity: dispenser: 100 pcs carton: 10 dispenser palet: 63 cartons Size: dispenser: L-24 cm, W-12cm, H-7 cm carton: L-36cm, W-25cm, H-25cm palet: L-120 cm, W-80 cm, H-190 cm
Manufacturer:	 Mercator Medical S.A. ul. H. Modrzejewskiej 30 31-327 Kraków, Poland